

• HLAscan v2.0

1. Updates

- (1) HLAscan v2.0 is compatible with human reference (hg19 (v37) and v38).
- (2) Bam file and Fastq files can get to as an HLAscan v2.0 input.
- (3) Increase the ease of use.
- (4) Optimizations

HLAscan v2.0 uses Suffix Array (SA) method instead of BWA-mem for alignment with IMGT/HLA Database.

- Time efficiency: about 10 times faster than previous version (ex: under a minute for HLA-typing (<100dp))
- Space efficiency: < 20MB disk (no create intermediate files)

2. Contents

- (1) DB
HLA-ALL.IMGT : IMGT/HLA database
- (2) Tool
hla_scan_r_v2.0

3. Usage

- (1) Download the DB & Tool
- (2) Read Manual
./hla_scan
- (3) Execute the tool
 - Bam Input
 - ./hla_scan -b [Bam file] -v [37/38] -d [DB] -g [Gene Symbol] -t [# of threads]
 - Fastq Input
 - ./hla_scan -l [fastq file] (-r [fastq file2]) -v [37/38] -d [DB] -g [Gene Symbol] -t [# of threads]

4. Commands

- (1) HLA typing (input: bam file)
./hla_scan_v5.0 -b [bam] -d [IMGT/HLA DB] -v [version (37 or 38)]

(2) HLA typing (input: fastq files)

```
./hla_scan_v5.0 -l [fastq] (-r [fastq2]) -d [IMGT/HLA DB]
```

5. Options

(1) Gene List -g [string], default=HLA-A

HLA-A, HLA-B, HLA-C, HLA-E, HLA-F, HLA-G, MICA, MICB,
HLA-DMA, HLA-DMB, HLA-DOA, HLA-DOB, HLA-DPA1, HLA-DPB1, HLA-DQA1,
HLA-DQB1, HLA-DRA, HLA-DRB1, HLA-DRB5, TAP1, TAP2

(2) Score Cutoff -s [int], default=50

(3) Constant using ScoreFunc -c [int], default=30

(4) # of threads -t [int], default=32